

## PS-007 Splendor Supreme


| PRODUCT SPECIFICATIONS                   |                                |
|------------------------------------------|--------------------------------|
| <b>Model :</b> Splendor Supreme          | <b>Rating :</b> N/A            |
| <b>Size :</b> HQ, LS, KS, DB, LD, QN, KN | <b>Category :</b> Adjustable   |
| <b>Mobility/Medical :</b> Yes            | <b>NDIS Support List :</b> Yes |

| PRODUCT SUMMARY                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Mattress Assembly:</b> Designed and manufactured in accordance with Pan Pacific Clinical Practice Guidelines for the Prevention and Management of Pressure Injury;</p> <ul style="list-style-type: none"> <li>Types of support surface (table 6.1) <ul style="list-style-type: none"> <li>Complies as a Reactive Combination using hybrid format comprising of High specification foam, Gel and Air.</li> </ul> </li> <li>Exceeds the minimum recommendations for high specification foam mattress (table 6.2) <ul style="list-style-type: none"> <li><b>Classification AS 2281-1993 of H, HR or LR = Yes</b></li> <li><b>Multilayer = Yes</b></li> <li><b>Depth = &gt; 200mm</b></li> <li><b>Mattress Assembly MVTR = Customisable</b> to end user requirements with fabric or vinyl cover</li> </ul> </li> <li>Designed to provide enhanced pressure care support for individuals with a Waterlow rating of 15-20 (High Risk)</li> </ul> <p>This mattress is primarily designed for mobility and medical purposes and is compatible with adjustable bed bases.</p> |



| PRODUCT TESTING                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Fabric Safety Standard :</b> Compliant to OEKO-TEX Standard 100, Product Class 1, SHKO 064515</p> <p><b>Foam Safety Standard :</b> Compliant to OEKO-TEX Standard 100, Product Class 1, SHAO 116007</p> <p><b>Clinical Pressure Mapping :</b> Tested in accordance with AP-006 Pressure Mapping Procedure</p> <p><b>Bamboo Fibre :</b> Certified EU Ecolabel for textile products AT/016/001</p> <p><b>Pressure Memory Foam :</b> Tested in accordance with Australian Standard AS 2282</p> |

| CLINICAL PRESSURE MAPPING              |                                     |
|----------------------------------------|-------------------------------------|
| <b>Model :</b> Splendor Supreme        | <b>Location:</b> Brisbane Australia |
| <b>Test Equipment :</b> XSensor Pro V8 | <b>Test Date :</b> 31/07/25         |
| <b>Bed Surface :</b> SmartFlex 3       | <b>Ambient Temperature :</b> 22°C   |
| <b>Units :</b> mmHg                    | <b>Tested By :</b> Trent Collins    |
| <b>Test Subject Height :</b> 160cm     | <b>Test Subject Weight :</b> 65kg   |



| XSENSOR PRO v8 DETAILS                  |                            |
|-----------------------------------------|----------------------------|
| <b>Sensor:</b> PX100:64. 160.02         | <b>Rows :</b> 160          |
| <b>Resolution :</b> 0.5                 | <b>Columns :</b> 64        |
| <b>Calibration :</b> ISO/IEC 17025:2017 | <b>Frame Rate:</b> 443 fps |

Avante uses the XSENSOR pressure mat for accurate measurement and image quality to analyse pressure data. Pressure Mapping quantifies the interface pressure distribution between the body and the mattress to assist in ensuring peak and average pressures are assessed to reduce the risk of pressure injuries.